

Product Datasheet

FABP4 (Phospho-Tyr20) antibody (orb156815)

Description	FABP4 (Phospho-Tyr20) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthesised phosphopeptide derived from human FABP4 around the phosphorylation site of Tyr20 (DD(p-Y)MK)
Target	phospho-FABP4 (Tyr20)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14.6 kDa
Uniprot ID	P15090
Expiration Date	12 months from date of receipt.