


Product Datasheet

Factor I heavy chain antibody (orb156819)

Description	Factor I heavy chain antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Complement factor I heavy chain (261-360/582 aa)
Target	Factor I heavy chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35/63 kDa
Uniprot ID	P05156
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa (human) ce...