



Product Datasheet

FLNB antibody (orb156897)

Description	FLNB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human FLNB/Filamin 3 (1551-1650/2602 aa)
Target	FLNB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	278 kDa
Uniprot ID	075369
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa Human Cell...



Western blot analysis of U251 Human Cell...