

Product Datasheet

FLNC (2233) antibody (orb156898)

Description	FLNC (2233) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human FLNC around the phosphorylation site of Ser2233 (LG(p-S)FG)
Target	phospho-FLNC (Ser2233)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	300 kDa
Uniprot ID	Q14315
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.