


Product Datasheet

Fos B (Phospho-Ser27) antibody (orb156921)

Description	Fos B (Phospho-Ser27) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Fos B around the phosphorylation site of Ser27 (QYL(p-S)SV)
Target	phospho-Fos B (Ser27)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P53539
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Hela cell
us...