

Product Datasheet

FOXB1 + FOXB2 antibody (orb156924)

Description	FOXB1 + FOXB2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human FOXB1 + FOXB2 (51-100/325 aa)
Target	FOXB1 + FOXB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35/45 kDa
Uniprot ID	Q5VYV0
Expiration Date	12 months from date of receipt.