

Product Datasheet

FOXD3 antibody (orb156928)

Description	FOXD3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human FOXD3 (151-230/478 aa)
Target	FOXD3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	P19099
Dilution Range	ICC=1:25
Expiration Date	12 months from date of receipt.