


Product Datasheet

FOXH1 antibody (orb156933)

Description	FOXH1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human FOXH1 (41-140/365 aa)
Target	FOXH1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	075593
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Bone
(Mouse);Primary:Anti-
F...