



Product Datasheet

FOXO4 (Phospho-Thr451) antibody (orb156939)

Description	FOXO4 (Phospho-Thr451) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human FOXO4 around the phosphorylation site of Thr451 (LG(p-T)PV)
Target	phospho-FOXO4 (Thr451)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	P98177
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Placenta ...



Western blot analysis of K562 Cell Lysat...