


Product Datasheet

Frizzled 8 antibody (orb156954)

Description	Frizzled 8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Frizzled 8 (61-160/694 aa)
Target	Frizzled 8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03%, Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	Q9H461
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human U87MG cell; Lane 2: Actin