



Product Datasheet

GABA Transporter 3 antibody (orb156998)

Description	GABA Transporter 3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human GABA Transporter 3/GAT 3 (531-632/632 aa)
Target	GABA Transporter 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71 kDa
Uniprot ID	P48066
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/test, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from s...



Western blot analysis of A431 human cell...