



Product Datasheet

GABRB2 antibody (orb157008)


Description	GABRB2 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GABRB2/GABA A Receptor beta 2 (251-350/512 aa)
Target	GABRB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P47870
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of
Mouse brain
using GABRB2...



IHC-P analysis of
rat brain using
GABRB2...



Western
blot: Lane
1: Mouse
Cerebellum; Lan...