



Product Datasheet

CD100 Antibody (FITC) (orb1570908)


Description	CD100, also known as semaphorin 4D, is a homodimerizing type I transmembrane glycoprotein containing...
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	PHA stimulated human PBL
Target	CD100
Preservatives	Stabilizing Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgM
Clonality	Monoclonal
Uniprot ID	Q92854
Entrez	10507
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD100
positive
CD56 ...



Flow
cytometry
surface
staining
pattern ...