

Product Datasheet

Mouse CD39L1 / ENTPD2 Protein, His Tag (active enzyme, MALS verified) [Out of stock] (orb1570926)

Description

CD39L1, also known as ENTPD2

ectonucleotidase belonging to the CD39 family. It i...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CD39L1/ENTPD2

Preservatives

20mM Tris,150mM NaCl,pH7.5

Form/Appearance

Liquid

Storage

-70°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.2 kDa. The protein migrates as 44-46 kDa and 64-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse CD39L1 / ENTPD2 Protein, His Tag (orb1570926) is expressed from CHO cells. It contains AA Cys 26 - Ser 462 (Accession # O55026-1).

MW

50.2 kDa

Uniprot ID
O55026
Expiration Date

6 months from date of receipt.