



Product Datasheet

PDGFR beta Antibody (orb1570943)


Description	PDGFR beta Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	residues 32-530 (ECD)
Target	PDGFR beta
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P05622
NCBI	NP_032835.2
Dilution Range	WB - 1:1,000; IP - 6 µl/1 mg lysate; ICC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections.
Expiration Date	12 months from date of receipt.



Detection of mouse PDGFR beta by western...



Detection of mouse PDGFR beta (shaded) i...



Detection of mouse PDGFR beta by western...