



Product Datasheet

GFPT1 antibody (orb157116)


Description	GFPT1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GFPT1 (601-699/699 aa)
Target	GFPT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79 kDa
Uniprot ID	Q06210
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Western blot analysis of hela lysate usi...



Western blot analysis of 293T cell lysat...