


Product Datasheet

GGT2 heavy chain antibody (orb157127)

Description	GGT2 heavy chain antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Gamma-glutamyltranspeptidase 2 heavy chain (101-200/569 aa)
Target	GGT2 heavy chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	P36268
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
293FT Cell
u...