

Product Datasheet

GGT7 light chain antibody (orb157133)

Description	GGT7 light chain antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human GGT7 light chain (581-662/662 aa)
Target	GGT7 light chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q9UJ14
Expiration Date	12 months from date of receipt.