



Product Datasheet

GJA8 antibody (orb157149)

Description	GJA8 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human GJA8/Connexin 50 (151-250/433 aa)
Target	GJA8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	P48165
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...



Immunohistochemical staining of Mouse em...