

Product Datasheet

Glucose 6 Phosphate Dehydrogenase antibody (orb157179)

Description	Glucose 6 Phosphate Dehydrogenase antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Glucose 6 Phosphate Dehydrogenase (351-450/515 aa)
Target	Glucose 6 Phosphate Dehydrogenase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	P11413
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.