

Product Datasheet

Golgin subfamily A member 4 antibody (orb157249)

Description	Golgin subfamily A member 4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Golgin subfamily A member 4/Golgin 245 (2061-2160/2230 aa)
Target	Golgin subfamily A member 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	261 kDa
Uniprot ID	Q13439
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.