

Product Datasheet

Human PD-1 Antibody (PE) (orb1572507)

Description	Human PD-1 Antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant Human PD-1 Protein
Target	Human PD-1
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q15116
Expiration Date	12 months from date of receipt.