



Product Datasheet

LGR7 antibody (orb157275)


Description	LGR7 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LGR7. Please contact us for the exact immunogen sequence. The peptide is available as orb12941.
Target	LGR7
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	57 kDa
Uniprot ID	Q6R6I6, Q6R6I7, Q9HBX9
NCBI	001253728, 001240657
Entrez	59350
Dilution Range	IF/ICC: 1:50-500, IHC-P: 1:50-500, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat ovary (lane 1)



human thyroid cancer (lane 5), 293T cell...



IHC-P image of pig spleen tissue using a...