



Product Datasheet

GPR1 antibody (orb157288)


Description	GPR1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human GPR1 (151-250/355 aa)
Target	GPR1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P46091
Expiration Date	12 months from date of receipt.



WB analysis of rat kidney (lane 1), rat ...



Immunohistochemical staining of paraffin...



IHC-P staining of rat kidney tissue usin...