


Product Datasheet

GPR125 antibody (orb157315)

Description	GPR125 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GPR125 (351-450/1321 aa)
Target	GPR125
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	142 kDa
Uniprot ID	Q8IWK6
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human br...