



Product Datasheet

GPR147 antibody (orb157323)


Description	GPR147 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GPR147/NPFF1 Receptor (151-260/430 aa)
Target	NPFF1 Receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q9GZQ6
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Sample:
U251(Human) Cell
Lysate at 30 ug...



Paraformaldehyde-
fixed, paraffin
embedde...



Paraformaldehyde-
fixed, paraffin
embedde...