



Product Datasheet

GPR177 antibody (orb157332)


Description	GPR177 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GPR177 (351-450/541 aa)
Target	GPR177
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62 kDa
Uniprot ID	Q5T9L3
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Brain lys...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat brai...