


Product Datasheet

GRAF1 antibody (orb157379)

Description	GRAF1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GRAF1 (1-100/814 aa)
Target	GRAF
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	92 kDa
Uniprot ID	Q9UNA1
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Cerebrum(Rat); Cerebrum(Mouse); Cerebrum(Human); Liver(Mouse); Liver(Human); Liver(Rat); Kidney(Mouse); Kidney(Human); Kidney(Rat); Control