

Product Datasheet

H3N2 Matrix Protein 2 antibody (orb157415)

Description	H3N2 Matrix Protein 2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human H3N2 Matrix Protein 2 (2-50/97 aa)
Target	H3N2 Matrix Protein 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P03452
Expiration Date	12 months from date of receipt.