

Product Datasheet

Histone H4 (Acetyl K5) Antibody (Cy5.5) (orb1574505)

Description	Histone H4 (Acetyl K5) Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PerCP/Cy5.5
Tested Applications	ICC
Immunogen	KLH conjugated synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Acetyl Lys5
Target	Histone H4 (Acetyl K5)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Expiration Date	12 months from date of receipt.