



Product Datasheet

HDAC4 (Phospho-Ser246) antibody (orb157455)

Description	HDAC4 (Phospho-Ser246) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human HDAC4 around the phosphorylation site of Ser246 (TA(p-S)EP)
Target	Phospho-HDAC4 (Ser246)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119 kDa
Uniprot ID	P56524
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Sample:
Hela(Human) Cell
Lysate at 30 ug...



Paraformaldehyde-
fixed, paraffin
embedde...