

Product Datasheet

HDAC4 (Phospho-Ser632) antibody (orb157456)

Description	HDAC4 (Phospho-Ser632) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human HDAC4 around the phosphorylation site of Ser632 (AQ(p-S)SP)
Target	phospho-HDAC4 (Ser632)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119 kDa
Uniprot ID	P56524
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.