

Product Datasheet

HDC2 antibody (orb157462)

Description	HDC2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human HDC2/PHC2 (651-750/858 aa)
Target	PHC2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	91 kDa
Uniprot ID	Q8IXK0
Expiration Date	12 months from date of receipt.