



Product Datasheet

HELZ antibody (orb157482)

Description	HELZ antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human HELZ (1-100/1942 aa)
Target	HELZ
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	219 kDa
Uniprot ID	P42694
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, Flow-Cyt=1µg/Test, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Flow cytometric analysis of Rsc96 cell U...