

Product Datasheet

CD90/Thy-1 Antibody (orb1574985)

Description	CD90/Thy-1 Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Human Thy-1 protein, C-His (HEK293) (20-141/161 aa)
Target	CD90/Thy-1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	12 kDa
Uniprot ID	P04216
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.