

Product Datasheet

CD11a Antibody (Cy3) (orb1574997)

Description	CD11a Antibody (Cy3)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy3
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human CD11a
Target	CD11a
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P20701
Dilution Range	Flow-Cyt=1ug/Test,ICC=1:50
Expiration Date	12 months from date of receipt.