

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid Antibody (AP) (orb1575089)

Description

SARS-CoV-2 (2019-nCoV)

Species/Host	Mouse
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.