



Product Datasheet

HEXB chain B antibody (orb157517)

Description	HEXB chain B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human HEXB (201-300/556 aa)
Target	HEXB chain B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22/50 kDa
Uniprot ID	P07686
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:HeLa(Human) Cell;Lung (Mous...



Western blot:Siha (Human) Cell;Bone (Mou...