



Product Datasheet

HGFA Inhibitor 2 antibody (orb157523)

Description	HGFA Inhibitor 2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human SPINT2 (31-130/252 aa)
Target	HGFA Inhibitor 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	O43291
Dilution Range	Flow-Cyt=1µg/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 Cell Human...



Siha Cell Human Lysate using HGFA Inhibi...