

Product Datasheet

HIV2 gp120 + gp160 antibody (orb157553)

Description	HIV2 gp120 + gp160 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from HIV2 gp120 + gp160 (311-380/858 aa)
Target	HIV2 gp120 + gp160
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q79670
Expiration Date	12 months from date of receipt.