



Product Datasheet

IDO antibody (orb157614)


Description	IDO antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IDO. Please contact us for the exact immunogen sequence. The peptide is available as orb12942.
Target	IDO
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for one month and at -20°C for a longer time. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	45 kDa
Uniprot ID	P14902, P28776, Q9ERD9
Entrez	3620
Dilution Range	IF/ICC: 1:50-400, IHC-P: 1:50-300, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



IF analysis of rat lung tissue using ant...



Immunohistochemical staining of paraffin...



IHC-P image of rat lung tissue using ant...