

Product Datasheet

Iduronate 2 sulfatase antibody (orb157616)

Description

Iduronate 2 sulfatase antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Iduronate 2 sulfatase. Please contact us for the exact immunogen sequence. The peptide is available as orb374816.

Target

Iduronate 2 sulfatase

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

62 kDa

Uniprot ID

[P22304](#)

Entrez



[3423](#)

Dilution Range


WB: 1:100-500, IHC-P: 1:100-500 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P image of rat pancreas tissue using...



IHC-P staining of mouse testis tissue us...



Immunohistochemical staining of rat kidn...