

Product Datasheet

phospho-VEGF receptor 2 (Tyr1059) Antibody (Cy5.5) (orb1576181)

Description	phospho-VEGF receptor 2 (Tyr1059) Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human VEGF receptor 2 around the phosphorylation site of Tyr1059
Target	phospho-VEGF receptor 2 (Tyr1059)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35968
Expiration Date	12 months from date of receipt.