

Product Datasheet

IFNAR1 (Phospho-Tyr466) antibody (orb157643)

Description	IFNAR1 (Phospho-Tyr466) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IFNAR1 around the phosphorylation site of Tyr466 (IN(p-Y)VF)
Target	phospho-IFNAR1 (Tyr466)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	P17181
Expiration Date	12 months from date of receipt.