

Product Datasheet

Netrin G1 ligand Antibody (BF488) (orb1576711)

Description	Netrin G1 ligand Antibody (BF488)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF488
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Netrin G1 ligand
Target	Netrin G1 ligand
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCJ2
Expiration Date	12 months from date of receipt.