



Product Datasheet

IL-19 antibody (orb157681)

Description	IL-19 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human IL-19/NG.1 (81-177/177 aa)
Target	IL-19
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	Q9UHD0
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Mouse Spleen; Lane
2: ...



Western blot: Jurkat
cell(human); Primary: ...