

Product Datasheet

IL2 Receptor beta (Phospho-Tyr364) antibody (orb157684)

Description	IL2 Receptor beta (Phospho-Tyr364) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IL2 Receptor beta around the phosphorylation site of Tyr364 (QG(p-Y)F)
Target	phospho-IL2 Receptor beta (Tyr364)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P14784
Expiration Date	12 months from date of receipt.