


Product Datasheet

KCNG3 antibody (orb157733)

Description	KCNG3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human KCNG3 (251-350/436 aa)
Target	KCNG3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q8TAE7
Dilution Range	WB=1:500-2000, IF=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human Jurkat cell; La...