


Product Datasheet

KMT6 (Phospho-Ser21) antibody (orb157761)

Description	KMT6 (Phospho-Ser21) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human KMT6/EZH2 around the phosphorylation site of Ser21 (VK(p-S)EY)
Target	phospho-KMT6 (Ser21)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	Q15910
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Paraformaldehyde-fixed, paraffin embedde...