

Product Datasheet

LPAAT-epsilon antibody (orb157812)

Description	LPAAT-epsilon antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human LPAAT-epsilon/AGPAT5 (188-285/364 aa)
Target	AGPAT5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	Q9NUQ2
Expiration Date	12 months from date of receipt.